

Antibodies Product Datasheet

Product Name:	FRMD3 Antibody
Catalog Number:	C43831
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FRMD3
Concentration:	1 mg/ml
Gene Name:	FRMD3 EPB41L4O
Protein Name:	FERM domain-containing protein 3 (Band 4.1-like protein 4O) (Ovary type protein 4.1) (4.1O)
Other Name:	
Human Swiss Prot No.:	A2A2Y4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FRMD3
Storage and Stability:	-20°C/1 year