

Antibodies Product Datasheet

Product Name:	FRMD5 Antibody
Catalog Number:	C43832
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FRMD5
Concentration:	1 mg/ml
Gene Name:	FRMD5
Protein Name:	FRMD5
Other Name:	FERM domain-containing protein 5
Human Swiss Prot No.:	Q7Z6J6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.182% sodium azide.
Specificity:	This antibody detects endogenous levels of FRMD5 at Human/Mouse
Storage and Stability:	-20°C/1 year