

Antibodies Product Datasheet

Product Name:	FTSJD1 Antibody
Catalog Number:	C43840
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FTSJ1
Concentration:	1 mg/ml
Gene Name:	FTSJD1 AFT
Protein Name:	FTSJ1
Other Name:	FtsJ methyltransferase domain-containing protein 1 (EC 2.1.1.-) (Protein adrift homolog)
Human Swiss Prot No.:	Q8IYT2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.334% sodium azide.
Specificity:	This antibody detects endogenous levels of FTSJ1 at Human/Mouse
Storage and Stability:	-20°C/1 year