

Antibodies Product Datasheet

Product Name:	FXL20 Antibody
Catalog Number:	C43847
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FXL20
Concentration:	1 mg/ml
Gene Name:	FBXL20 FBL2
Protein Name:	F-box/LRR-repeat protein 20 (F-box and leucine-rich repeat protein 20) (F-box/LRR-repeat protein 2-like)
Other Name:	
Human Swiss Prot No.:	Q96IG2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FXL20
Storage and Stability:	-20°C/1 year