

Antibodies Product Datasheet

Product Name:	FXYD4 Antibody
Catalog Number:	C43848
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FXYD4
Concentration:	1 mg/ml
Gene Name:	FXYD4 UNQ526/PRO1069
Protein Name:	FXYD4
Other Name:	
Human Swiss Prot No.:	P59646
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FXYD4 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year