

Antibodies Product Datasheet

Product Name:	G3PT Antibody
Catalog Number:	C43851
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	GAPDHS GAPD2 GAPDH2 GAPDS HSD-35 HSD35
Protein Name:	Glyceraldehyde-3-phosphate dehydrogenase, testis-specific (EC 1.2.1.12) (Spermatogenic cell-specific glyceraldehyde 3-phosphate dehydrogenase 2) (GAPDH-2) (Spermatogenic glyceraldehyde-3-phosphate dehydrogenase)
Other Name:	
Human Swiss Prot No.:	O14556
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	G3PT Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year