

Antibodies Product Datasheet

Product Name:	G3ST1 Antibody
Catalog Number:	C43852
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human G3ST1
Concentration:	1 mg/ml
Gene Name:	GAL3ST1 CST
Protein Name:	Galactosylceramide sulfotransferase (GalCer sulfotransferase) (EC 2.8.2.11) (3'-phosphoadenosine-5'-phosphosulfate:GalCer sulfotransferase) (3'-phosphoadenylylsulfate:galactosylceramide 3'-sulfotransferase) (Cerebroside sulfotransferase)
Other Name:	
Human Swiss Prot No.:	Q99999
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human G3ST1
Storage and Stability:	-20°C/1 year