

Antibodies Product Datasheet

Product Name:	GALT2 Antibody
Catalog Number:	C43860
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GALT2
Concentration:	1 mg/ml
Gene Name:	GALNT2
Protein Name:	GALT2
Other Name:	Polypeptide N-acetylgalactosaminyltransferase 2 (EC 2.4.1.41) (Polypeptide GalNAc transferase 2) (GalNAc-T2) (pp-GaNTase 2) (Protein-UDP acetylgalactosaminyltransferase 2) (UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 2) [Cleaved into: Polypep
Human Swiss Prot No.:	Q10471
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.194% sodium azide.
Specificity:	This antibody detects endogenous levels of GALT2 at Human/Mouse
Storage and Stability:	-20°C/1 year