

Antibodies Product Datasheet

Product Name:	GBLP Antibody
Catalog Number:	C43873
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	GNB2L1 HLC7 PIG21
Protein Name:	Guanine nucleotide-binding protein subunit beta-2-like 1 (Cell proliferation-inducing gene 21 protein) (Guanine nucleotide-binding protein subunit beta-like protein 12.3) (Human lung cancer oncogene 7 protein) (HLC-7) (Receptor for activated C kinase) (Re
Other Name:	
Human Swiss Prot No.:	P63244
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GBLP Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year