

Antibodies Product Datasheet

Product Name:	GCNT4 Antibody
Catalog Number:	C43888
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	GCNT4
Protein Name:	Beta-1,3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 4 (EC 2.4.1.102) (Core 2-branching enzyme 3) (Core2-GlcNAc-transferase 3) (C2GnT3)
Other Name:	
Human Swiss Prot No.:	Q9P109
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GCNT4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year