

Antibodies Product Datasheet

Product Name:	GDN Antibody
Catalog Number:	C43891
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	SERPINE2 PI7 PN1
Protein Name:	Glia-derived nexin (GDN) (Peptidase inhibitor 7) (PI-7) (Protease nexin 1) (PN-1) (Protease nexin I) (Serpine E2)
Other Name:	
Human Swiss Prot No.:	P07093
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GDN Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year