

Antibodies Product Datasheet

Product Name:	GDS1 Antibody
Catalog Number:	C43893
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 230-310
Concentration:	1 mg/ml
Gene Name:	RAP1GDS1
Protein Name:	Rap1 GTPase-GDP dissociation stimulator 1 (Exchange factor smgGDS) (SMG GDS protein) (SMG P21 stimulatory GDP/GTP exchange protein)
Other Name:	
Human Swiss Prot No.:	P52306
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GDS1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year