

Antibodies Product Datasheet

Product Name:	GEMI2 Antibody
Catalog Number:	C43896
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	GEMIN2 SIP1
Protein Name:	Gem-associated protein 2 (Gemin-2) (Component of gems 2) (Survival of motor neuron protein-interacting protein 1) (SMN-interacting protein 1)
Other Name:	
Human Swiss Prot No.:	O14893
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GEMI2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year