

Antibodies Product Datasheet

Product Name:	GEMI4 Antibody
Catalog Number:	C43897
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GEMI4
Concentration:	1 mg/ml
Gene Name:	GEMIN4
Protein Name:	Gem-associated protein 4 (Gemin-4) (Component of gems 4) (p97)
Other Name:	
Human Swiss Prot No.:	P57678
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GEMI4
Storage and Stability:	-20°C/1 year