

Antibodies Product Datasheet

Product Name:	GFRA2 Antibody
Catalog Number:	C43903
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 320-400
Concentration:	1 mg/ml
Gene Name:	GFRA2 GDNFRB RETL2 TRNR2
Protein Name:	GNDF family receptor alpha-2 (GNDF receptor alpha-2) (GDNFR-alpha-2) (GFR-alpha-2) (GNDF receptor beta) (GDNFR-beta) (Neurturin receptor alpha) (NRTNR-alpha) (NTNR-alpha) (RET ligand 2) (TGF-beta-related neurotrophic factor receptor 2)
Other Name:	
Human Swiss Prot No.:	O00451
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GFRA2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year