

Antibodies Product Datasheet

Product Name:	GGT3 Antibody
Catalog Number:	C43911
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GGT3
Concentration:	1 mg/ml
Gene Name:	GGT3P GGT3
Protein Name:	GGT3
Other Name:	Putative gamma-glutamyltranspeptidase 3 (GGT 3) (EC 2.3.2.2) (Gamma-glutamyltransferase 3) (Glutathione hydrolase 3) (EC 3.4.19.13) [Cleaved into: Putative gamma-glutamyltranspeptidase 3 heavy chain; Putative gamma-glutamyltranspeptidase 3 light chain]
Human Swiss Prot No.:	A6NGU5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.154% sodium azide.
Specificity:	This antibody detects endogenous levels of GGT3 at Human
Storage and Stability:	-20°C/1 year