

Antibodies Product Datasheet

Product Name:	GLP2R Antibody
Catalog Number:	C43920
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:500-1000, IHC 1:50-100
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 180-260
Concentration:	1 mg/ml
Gene Name:	GLP2R
Protein Name:	Glucagon-like peptide 2 receptor (GLP-2 receptor) (GLP-2-R) (GLP-2R)
Other Name:	
Human Swiss Prot No.:	O95838
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GLP2R Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year