

Antibodies Product Datasheet

Product Name:	GP183 Antibody
Catalog Number:	C43963
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 110-190
Concentration:	1 mg/ml
Gene Name:	GPR183 EBI2
Protein Name:	G-protein coupled receptor 183 (Epstein-Barr virus-induced G-protein coupled receptor 2) (EBI2) (EBV-induced G-protein coupled receptor 2)
Other Name:	
Human Swiss Prot No.:	P32249
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GP183 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year