

Antibodies Product Datasheet

Product Name:	GPER Antibody
Catalog Number:	C43968
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 300-380
Concentration:	1 mg/ml
Gene Name:	GPER CEPR CMKRL2 DRY12 GPR30
Protein Name:	G-protein coupled estrogen receptor 1 (Chemoattractant receptor-like 2) (Flow-induced endothelial G-protein coupled receptor 1) (FEG-1) (G-protein coupled receptor 30) (GPCR-BR) (IL8-related receptor DRY12) (Lymphocyte-derived G-protein coupled receptor)
Other Name:	
Human Swiss Prot No.:	Q99527
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GPER Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year