

Antibodies Product Datasheet

Product Name:	GPR98 Antibody
Catalog Number:	C43989
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	GPR98 KIAA0686 KIAA1943 MASS1 VLGR1
Protein Name:	G-protein coupled receptor 98 (Monogenic audiogenic seizure susceptibility protein 1 homolog) (Usher syndrome type-2C protein) (Very large G-protein coupled receptor 1)
Other Name:	
Human Swiss Prot No.:	Q8WVG9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GPR98 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year