

Antibodies Product Datasheet

Product Name:	GPSM3 Antibody
Catalog Number:	C43991
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	GPSM3 AGS4 C6orf9 G18
Protein Name:	G-protein-signaling modulator 3 (Activator of G-protein signaling 4) (G18.1b) (Protein G18)
Other Name:	
Human Swiss Prot No.:	Q9Y4H4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GPSM3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year