

Antibodies Product Datasheet

Product Name:	GPX3 Antibody
Catalog Number:	C43996
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	GPX3 GPXP
Protein Name:	Glutathione peroxidase 3 (GPx-3) (GSHPx-3) (EC 1.11.1.9) (Extracellular glutathione peroxidase) (Plasma glutathione peroxidase) (GPx-P) (GSHPx-P)
Other Name:	
Human Swiss Prot No.:	P22352
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GPX3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year