

Antibodies Product Datasheet

Product Name:	GRAM2 Antibody
Catalog Number:	C43999
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GRAM2
Concentration:	1 mg/ml
Gene Name:	GRAMD2
Protein Name:	GRAM2
Other Name:	GRAM domain-containing protein 2
Human Swiss Prot No.:	Q8IU3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.54% sodium azide.
Specificity:	This antibody detects endogenous levels of GRAM2 at Human/Mouse
Storage and Stability:	-20°C/1 year