

Antibodies Product Datasheet

Product Name:	GRP2 Antibody
Catalog Number:	C44008
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GRP2
Concentration:	1 mg/ml
Gene Name:	RASGRP2 CDC25L MCG7
Protein Name:	RAS guanyl-releasing protein 2 (Calcium and DAG-regulated guanine nucleotide exchange factor I) (CaDAG-GEFI) (Cdc25-like protein) (hCDC25L) (F25B3.3 kinase-like protein)
Other Name:	
Human Swiss Prot No.:	Q7LDG7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GRP2
Storage and Stability:	-20°C/1 year