

Antibodies Product Datasheet

Product Name:	GTD2A Antibody
Catalog Number:	C44017
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GTD2A
Concentration:	1 mg/ml
Gene Name:	GTF2IRD2 GTF2IRD2A
Protein Name:	GTD2A
Other Name:	General transcription factor II-I repeat domain-containing protein 2A (GTF2I repeat domain-containing protein 2A) (Transcription factor GTF2IRD2-alpha)
Human Swiss Prot No.:	Q86UP8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.202% sodium azide.
Specificity:	This antibody detects endogenous levels of GTD2A at Human
Storage and Stability:	-20°C/1 year