

Antibodies Product Datasheet

Product Name:	GTD2B Antibody
Catalog Number:	C44018
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GTD2B
Concentration:	1 mg/ml
Gene Name:	GTF2IRD2B
Protein Name:	GTD2B
Other Name:	General transcription factor II-I repeat domain-containing protein 2B (GTF2I repeat domain-containing protein 2B) (Transcription factor GTF2IRD2-beta)
Human Swiss Prot No.:	Q6EKJ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.203% sodium azide.
Specificity:	This antibody detects endogenous levels of GTD2B at Human
Storage and Stability:	-20°C/1 year