

Antibodies Product Datasheet

Product Name:	H2B1J Antibody
Catalog Number:	C44034
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human H2B1J
Concentration:	1 mg/ml
Gene Name:	HIST1H2BJ H2BFR
Protein Name:	Histone H2B type 1-J (Histone H2B.1) (Histone H2B.r) (H2B/r)
Other Name:	
Human Swiss Prot No.:	P06899
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human H2B1J
Storage and Stability:	-20°C/1 year