

Antibodies Product Datasheet

Product Name:	HBAT Antibody
Catalog Number:	C44050
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HBQ1
Protein Name:	Hemoglobin subunit theta-1 (Hemoglobin theta-1 chain) (Theta-1-globin)
Other Name:	
Human Swiss Prot No.:	P09105
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HBAT Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year