

Antibodies Product Datasheet

Product Name:	HCAR2 Antibody
Catalog Number:	C44058
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 260-340
Concentration:	1 mg/ml
Gene Name:	HCAR2 GPR109A HCA2 HM74A NIACR1
Protein Name:	Hydroxycarboxylic acid receptor 2 (G-protein coupled receptor 109A) (G-protein coupled receptor HM74A) (Niacin receptor 1) (Nicotinic acid receptor)
Other Name:	
Human Swiss Prot No.:	Q8TDS4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HCAR2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year