

Antibodies Product Datasheet

Product Name:	HDHD1 Antibody
Catalog Number:	C44063
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HDHD1
Concentration:	1 mg/ml
Gene Name:	HDHD1 DXF68S1E FAM16AX GS1 HDHD1A
Protein Name:	Pseudouridine-5'-monophosphatase (5'-PsiMPase) (EC 3.1.3.n6) (Haloacid dehalogenase-like hydrolase domain-containing protein 1) (Haloacid dehalogenase-like hydrolase domain-containing protein 1A) (Protein GS1)
Other Name:	
Human Swiss Prot No.:	Q08623
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HDHD1
Storage and Stability:	-20°C/1 year