

Antibodies Product Datasheet

Product Name:	HECW1 Antibody
Catalog Number:	C44067
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HECW1 KIAA0322 NEDL1
Protein Name:	E3 ubiquitin-protein ligase HECW1 (EC 6.3.2.-) (HECT, C2 and WW domain-containing protein 1) (NEDD4-like E3 ubiquitin-protein ligase 1) (hNEDL1)
Other Name:	
Human Swiss Prot No.:	Q76N89
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HECW1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year