

Antibodies Product Datasheet

Product Name:	HEN2 Antibody
Catalog Number:	C44074
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	NHLH2 BHLHA34 HEN2 KIAA0490
Protein Name:	Helix-loop-helix protein 2 (HEN-2) (Class A basic helix-loop-helix protein 34) (bHLHa34) (Nescient helix loop helix 2) (NSCL-2)
Other Name:	
Human Swiss Prot No.:	Q02577
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HEN2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year