

Antibodies Product Datasheet

Product Name:	HERC1 Antibody
Catalog Number:	C44076
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HERC1
Concentration:	1 mg/ml
Gene Name:	HERC1
Protein Name:	Probable E3 ubiquitin-protein ligase HERC1 (EC 6.3.2.-) (HECT domain and RCC1-like domain-containing protein 1) (p532) (p619)
Other Name:	
Human Swiss Prot No.:	Q15751
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HERC1
Storage and Stability:	-20°C/1 year