

Antibodies Product Datasheet

Product Name:	HES3 Antibody
Catalog Number:	C44079
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 20-100
Concentration:	1 mg/ml
Gene Name:	HES3 BHLHB43
Protein Name:	Transcription factor HES-3 (Class B basic helix-loop-helix protein 43) (bHLHb43) (Hairy and enhancer of split 3)
Other Name:	
Human Swiss Prot No.:	Q5TGS1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HES3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year