

Antibodies Product Datasheet

Product Name:	HEY1 Antibody
Catalog Number:	C44084
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HEY1 BHLHB31 CHF2 HERP2 HESR1 HRT1
Protein Name:	Hairy/enhancer-of-split related with YRPW motif protein 1 (Cardiovascular helix-loop-helix factor 2) (CHF-2) (Class B basic helix-loop-helix protein 31) (bHLHb31) (HES-related repressor protein 1) (Hairy and enhancer of split-related protein 1) (HESR-1) (
Other Name:	
Human Swiss Prot No.:	Q9Y5J3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HEY1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year