

Antibodies Product Datasheet

Product Name:	HIC2 Antibody
Catalog Number:	C44090
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HIC2 HRG22 KIAA1020 ZBTB30
Protein Name:	Hypermethylated in cancer 2 protein (Hic-2) (HIC1-related gene on chromosome 22 protein) (Hic-3) (Zinc finger and BTB domain-containing protein 30)
Other Name:	
Human Swiss Prot No.:	Q96JB3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HIC2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year