

Antibodies Product Datasheet

Product Name:	HIPK3 Antibody
Catalog Number:	C44096
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 770-850
Concentration:	1 mg/ml
Gene Name:	HIPK3 DYRK6 FIST3 PKY
Protein Name:	Homeodomain-interacting protein kinase 3 (EC 2.7.11.1) (Androgen receptor-interacting nuclear protein kinase) (ANPK) (Fas-interacting serine/threonine-protein kinase) (FIST) (Homolog of protein kinase YAK1)
Other Name:	
Human Swiss Prot No.:	Q9H422
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HIPK3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year