

Antibodies Product Datasheet

Product Name:	HLAG Antibody
Catalog Number:	C44099
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HLA-G HLA-6.0 HLAG
Protein Name:	HLA class I histocompatibility antigen, alpha chain G (HLA G antigen) (MHC class I antigen G)
Other Name:	
Human Swiss Prot No.:	P17693
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HLAG Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year