

Antibodies Product Datasheet

Product Name:	HPGDS Antibody
Catalog Number:	C44122
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HPGDS GSTS PGDS PTGDS2
Protein Name:	Hematopoietic prostaglandin D synthase (H-PGDS) (EC 5.3.99.2) (GST class-sigma) (Glutathione S-transferase) (EC 2.5.1.18) (Glutathione-dependent PGD synthase) (Glutathione-requiring prostaglandin D synthase) (Prostaglandin-H2 D-isomerase)
Other Name:	
Human Swiss Prot No.:	O60760
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HPGDS Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year