

Antibodies Product Datasheet

Product Name:	HRSL2 Antibody
Catalog Number:	C44130
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HRSL2
Concentration:	1 mg/ml
Gene Name:	HRASLS2
Protein Name:	HRAS-like suppressor 2 (EC 2.3.1.-) (EC 3.1.1.-)
Other Name:	
Human Swiss Prot No.:	Q9NWW9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HRSL2
Storage and Stability:	-20°C/1 year