

Antibodies Product Datasheet

Product Name:	HS3S1 Antibody
Catalog Number:	C44133
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HS3S1
Concentration:	1 mg/ml
Gene Name:	HS3ST1 3OST 3OST1
Protein Name:	Heparan sulfate glucosamine 3-O-sulfotransferase 1 (EC 2.8.2.23) (Heparan sulfate D-glucosaminyl 3-O-sulfotransferase 1) (3-OST-1) (Heparan sulfate 3-O-sulfotransferase 1) (h3-OST-1)
Other Name:	
Human Swiss Prot No.:	O14792
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HS3S1
Storage and Stability:	-20°C/1 year