

Antibodies Product Datasheet

Product Name:	HS74L Antibody
Catalog Number:	C44135
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 450-530
Concentration:	1 mg/ml
Gene Name:	HSPA4L APG1 OSP94
Protein Name:	Heat shock 70 kDa protein 4L (Heat shock 70-related protein APG-1) (Osmotic stress protein 94)
Other Name:	
Human Swiss Prot No.:	O95757
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HS74L Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year