

Antibodies Product Datasheet

Product Name:	HSH2D Antibody
Catalog Number:	C44140
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HSH2D
Concentration:	1 mg/ml
Gene Name:	HSH2D ALX
Protein Name:	HSH2D
Other Name:	Hematopoietic SH2 domain-containing protein (Hematopoietic SH2 protein) (Adaptor in lymphocytes of unknown function X)
Human Swiss Prot No.:	Q96JZ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.18% sodium azide.
Specificity:	This antibody detects endogenous levels of HSH2D at Human/Mouse
Storage and Stability:	-20°C/1 year