

Antibodies Product Datasheet

Product Name:	HXB8 Antibody
Catalog Number:	C44160
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	HOXB8 HOX2D
Protein Name:	Homeobox protein Hox-B8 (Homeobox protein Hox-2.4) (Homeobox protein Hox-2D)
Other Name:	
Human Swiss Prot No.:	P17481
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HXB8 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year