

Antibodies Product Datasheet

Product Name:	I10R2 Antibody
Catalog Number:	C44170
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 240-320
Concentration:	1 mg/ml
Gene Name:	IL10RB CRFB4 D21S58 D21S66
Protein Name:	Interleukin-10 receptor subunit beta (IL-10 receptor subunit beta) (IL-10R subunit beta) (IL-10RB) (Cytokine receptor class-II member 4) (Cytokine receptor family 2 member 4) (CRF2-4) (Interleukin-10 receptor subunit 2) (IL-10R subunit 2) (IL-10R2) (CD an
Other Name:	
Human Swiss Prot No.:	Q08334
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	I10R2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year