

Antibodies Product Datasheet

Product Name:	I17RD Antibody
Catalog Number:	C44172
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	IL17RD IL17RLM SEF UNQ6115/PRO20026
Protein Name:	Interleukin-17 receptor D (IL-17 receptor D) (IL-17RD) (IL17Rhom) (Interleukin-17 receptor-like protein) (Sef homolog) (hSef)
Other Name:	
Human Swiss Prot No.:	Q8NFM7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	I17RD Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year