

Antibodies Product Datasheet

Product Name:	I22R1 Antibody
Catalog Number:	C44174
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	IL22RA1 IL22R
Protein Name:	Interleukin-22 receptor subunit alpha-1 (IL-22 receptor subunit alpha-1) (IL-22R-alpha-1) (IL-22RA1) (Cytokine receptor class-II member 9) (Cytokine receptor family 2 member 9) (CRF2-9) (ZcytoR11)
Other Name:	
Human Swiss Prot No.:	Q8N6P7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	I22R1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year