

Antibodies Product Datasheet

Product Name:	I23O1 Antibody
Catalog Number:	C44175
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	IDO1 IDO INDO
Protein Name:	Indoleamine 2,3-dioxygenase 1 (IDO-1) (EC 1.13.11.52) (Indoleamine-pyrrole 2,3-dioxygenase)
Other Name:	
Human Swiss Prot No.:	P14902
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	I23O1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year