

Antibodies Product Datasheet

Product Name:	IBP5 Antibody
Catalog Number:	C44177
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 90-170
Concentration:	1 mg/ml
Gene Name:	IGFBP5 IBP5
Protein Name:	Insulin-like growth factor-binding protein 5 (IBP-5) (IGF-binding protein 5) (IGFBP-5)
Other Name:	
Human Swiss Prot No.:	P24593
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IBP5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year