

Antibodies Product Datasheet

Product Name:	IBPL1 Antibody
Catalog Number:	C44178
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	IGFBPL1 IGFBPRP4
Protein Name:	Insulin-like growth factor-binding protein-like 1 (IGFBP-related protein 10) (Insulin-like growth factor-binding-related protein 4) (IGFBP-rP4)
Other Name:	
Human Swiss Prot No.:	Q8WX77
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IBPL1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year